

Interactive
Comment

Interactive comment on “Technical Note: A fully automated purge and trap-GC-MS system for quantification of volatile organic compound (VOC) fluxes between the ocean and atmosphere” by S. J. Andrews et al.

S. J. Andrews et al.

stephen.andrews@york.ac.uk

Received and published: 21 February 2015

Thank you for your positive comments. On page 5 (2983), line 7 we state the chromatographic settings used predominantly in this work and to produce the chromatogram in Fig.2: “The gas sample is pre-concentrated onto the cooled (-30°C) absorbent trap (Tenax TA) of the TD unit and desorbed at 250°C onto the GC column (60 m 320 μm I.D. 1 μm 5% phenyl film, SGE). The column was held at 40°C for 2 min, ramped at 20°Cmin⁻¹ to 207°C then at 100°Cmin⁻¹ to 250°C.”

Full Screen / Esc

Printer-friendly Version

Interactive Discussion

Discussion Paper



Backflushing was not used on the water instrument during any of the cruises but has been tested in the lab.

Interactive comment on Ocean Sci. Discuss., 11, 2979, 2014.

OSD

11, C1350–C1351, 2015

Interactive
Comment

Full Screen / Esc

Printer-friendly Version

Interactive Discussion

Discussion Paper



C1351